



20107002C
Xiros Mikro
Freeze Dryer
User Manual

Table of Contents

I. Disclaimer & Warranty	4
LIMITED WARRANTY & DISCLAIMER	4
LIMITATION OF LIABILITY.....	6
II. Safety	7
2.1 Panel Warning	9
2.2 Lifting and Moving.....	9
III. Sécurité.....	10
3.1 Avertissements du panneau.....	12
3.2 Levage et déplacement.....	12
IV. Introduction	13
4.1 After Sales Support.....	13
4.2 Proper Use	13
4.3 Product Disposal.....	13
4.4 Warranty	13
4.5 Customer Responsibilities & Warranty Conditions	14
V. Technical Specifications.....	15
VI. Precautions.....	16
VII. Components Overview.....	17
VIII. Software Operation	19
8.1 AutoDry.....	20
8.1.1 AutoDry Process	21
8.1.2 Editing AutoDry Parameters	23
8.2 Program.....	26
8.2.1 AutoDry/Program Stages.....	28
8.2.2 Writing a Program & Parameters Explained	29
8.2.2.1 Freezing Stage	29
8.2.2.2 Evacuation Stage.....	30
8.2.2.3 Drying Stage.....	31
8.2.2.4 Storage stage	32
8.2.3 Renaming a Program	33
8.2.4 Loading a Preset Program.....	34
8.2.5 Starting a Program	36
8.2.6 Program Vacuum Oven (VO).....	40
8.3 Extra Drying	41
8.4 Power Recovery	43
8.5 Manual Control	44

8.5.1	Defrost Mode.....	46
8.6	System Settings.....	49
8.7	History Screen	51
8.7.1	History Export.....	52
IX.	Fault Diagnosis.....	53
X.	Cleaning and Maintenance	54
10.1	Cleaning	54
10.2	Maintenance.....	54

I. Disclaimer & Warranty

- The intellectual property rights for this manual belongs to Holland Green Science Corporation (hereinafter referred to as “the Company”).
- Products of the company are under the patent protection of China and other countries and regions (inclusive of obtained patents and patents currently under application).
- The company reserves the right to change product specifications without notice.
- The contents of this manual are subject to change without notice.
- The manual’s content must not be copied, extracted, or modified in any form without the Company’s permission.

LIMITED WARRANTY & DISCLAIMER

Products sold by the Company are intended solely for the use(s) specified by the manufacturer and are not to be used for any other purposes including, but not limited to, unauthorized commercial purposes. Buyer warrants and represents that Buyer will properly test, use, qualify, and/or validate each Product for intended use, and, to the extent authorized, manufacture and market any final articles made from the Company’s Products in accordance with the practices of a reasonable person who is an expert in the field, and in strict compliance with all applicable national, state and local laws and regulations. Buyer acknowledges and understands that the purchase of Products only conveys to Buyer the non-transferable right for only the Buyer to use the Products purchased in compliance with the applicable manufacturer’s Published User Statement, Limited User Statement, or Limited License, if any. While this manual contains certain representations of the potential dangers from the instrument or equipment, Buyer is solely responsible for conducting any research necessary to learn the hazards involved for any of its intended uses of the instrument or equipment purchased from the Company and to properly warn its customers, employees and other personnel who may be exposed to such instrument or equipment of any risks involved in using or handling the instrument or equipment.

The Company warrants that this product will operate or perform substantially in conformance with the Company’s published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended storage and usage by properly trained personnel, for the duration of the manufacturer’s warranty period set forth in this product documentation, published specifications or package inserts. If a manufacturer’s warranty period is not specified in the product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of Bill of Sale to Buyer for equipment, and the earlier of any stated expiration dates or ninety (90) days for all other products (the “Warranty Period”). The Company agrees during the Warranty Period, to repair or replace, at the Company’s option, defective Products so as to cause the same to operate in substantial conformance with said published specifications; provided that Buyer shall (a) promptly notify the Company in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after the Company’s review, the Company will provide Buyer with service data and/or a Return Material Authorization (“RMA”), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to the Company with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of the Company. All replaced parts shall become the property of the Company. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Terms

& Conditions agreed to by the purchase of this product. Consumables are expressly excluded from this warranty.

Notwithstanding the foregoing, Products supplied by the Company but not manufactured by the Company are not warranted by the Company, and the Company disclaims all warranties associated with the accuracy of any original manufacturer or third-party supplier Product information, but the Company agrees to assign to Buyer any warranty rights in such Product that the Company may have from the original manufacturer or third-party supplier, to the extent such assignment is allowed by such original manufacturer or third-party supplier.

In no event shall the Company have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Product in a manner for which they were not designed, (v) causes external to the Product such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Product or (vii) use of the Product in combination with equipment or software not supplied by the Company. If the Company determines that Product for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse the Company for all costs of investigating and responding to such request at the Company's then prevailing time and materials rates. If the Company provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay the Company therefor at the Company's then prevailing time and materials rates.

ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, RELOCATION OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN THE COMPANY WITHOUT THE COMPANY'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY THE COMPANY, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE PRODUCT.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, THE COMPANY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCTS, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR USE. THE COMPANY DOES NOT WARRANT THAT THE PRODUCTS ARE ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

NOTHING IN THIS MANUAL OR IN THE TERMS AND AGREEMENT AGREED TO AT THE TIME OF PURCHASE SHALL BE CONSTRUED TO PROVIDE FOR THE TRANSFER, ASSIGNMENT OR ALIENATION OF THE LIMITED WARRANTIES PROVIDED BY THE COMPANY HEREIN TO ANY OTHER PARTY OTHER THAN BUYER, INCLUDING BUT NOT LIMITED TO, SUCCESSORS, HEIRS, SUBSIDIARIES OR ANY OTHER THIRD PERSON.

THE WARRANTIES IDENTIFIED ABOVE ARE THE SOLE AND EXCLUSIVE WARRANTIES WITH RESPECT TO THIS PRODUCT AND ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. ALL OTHER WARRANTIES ARE EXPRESSLY DISCLAIMED, INCLUDING, WITHOUT LIMITATION, THOSE OF NON-INFRINGEMENT, MERCHANTABILITY, OR FITNESS FOR A SPECIFIC PURPOSE OR USE (WHETHER ARISING FROM STATUTE, OR OTHERWISE IN LAW, OR FROM A COURSE OF

DEALING, OR USAGE OF TRADE). THIS WARRANTY IS NOT TRANSFERABLE FROM THE ORIGINAL PURCHASER TO A SUBSEQUENT OWNER. FURTHER, THE COMPANY IS NOT LIABLE IN CASES OF DELIBERATE, NEGLIGENT OR ACCIDENTAL MISUSE OF THE PRODUCTS, USE WITH INAPPROPRIATE REAGENTS OR CONSUMABLES, DAMAGE CAUSED BY DISASTER, REPAIR OR MODIFICATIONS DONE BY ANYONE OTHER THAN THE COMPANY. THE COMPANY DOES NOT WARRANT THAT THE PRODUCT WILL NOT INFRINGE THE INTELLECTUAL PROPERTY RIGHTS OF A THIRD PARTY EITHER ALONE OR IN COMBINATION WITH OTHER PRODUCTS OR IN THE OPERATION OF ANY PROCESS. THE COMPANY'S TOTAL LIABILITY FOR BREACH OF THESE WARRANTIES SHALL BE LIMITED TO THE GROSS PURCHASE PRICE FOR THE PRODUCT AS STATED IN THE RELEVANT INVOICE OR ORDER UNDER WHICH THE ACTIVITY GIVING RISE TO LIABILITY ARISES. UNDER NO CIRCUMSTANCES SHALL THE COMPANY BE LIABLE FOR INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES (INCLUDING LOSS OF BUSINESS OR PROFITS) UNDER ANY THEORY OF LIABILITY HEREUNDER.

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

LIMITATION OF LIABILITY

NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF THE COMPANY UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF THE COMPANY FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER THIS SECTION) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THERETOFORE PAID BY BUYER TO THE COMPANY WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE HUNDRED THOUSAND DOLLARS (\$100,000), WHICHEVER IS LESS.

With respect to any software products incorporated in or forming a part of the Product hereunder Buyer intends and agrees that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "Buyer" or similar or derivative words are understood and agreed to mean "licensee". Notwithstanding anything to the contrary contained herein, The COMPANY or its licensor, as the case may be, retains all rights and interest in software products provided hereunder. See the Terms & Conditions of purchase for details regarding such license.

II. Safety



! DANGER (may cause serious damage to property and or casualties)

1. No user serviceable parts inside. Service by authorized personnel only.
2. Read this manual prior to operating the freeze dryer and observe and verify that all connections are properly secured.
3. Utility requirements: Ensure that the electrical connections meet the local standards and that the power supply is compatible with the freeze dryer nameplate. The power source must be properly grounded.
4. Do not operate the freeze dryer if there is any damage to the power cables.
5. The chamber and shelves may be extremely cold or hot –Check temperature before opening the chamber or touching the chamber 's internal components, to avoid injury from frostbite or burns.
6. Do not force open the chamber door while the chamber is under vacuum. Doing so may cause serious injury and cause significant damage to the freeze dryer.
7. Do not expose any body parts or entire bodies to vacuum.
8. The freeze dryer has been exclusively designed for freeze drying of non-flammable, non-corrosive, and non-hazardous solids or liquid containing products.
9. There is no user replaceable battery.

! WARNING (may cause property damage or personal injury)



1. The freeze dryer may not be operated prior to reading this manual.
2. The freeze dryer must only be opened by authorized service personnel.
3. Do not place heavy objects on top of or stand on the freeze dryer.
4. It is prohibited to install any unauthorized software onto the freeze dryer.
5. Keep hands away from the chamber door during pump down, to prevent crushing injuries.
6. Ensure all ventilation openings are not obstructed, and air can flow freely.
7. Do not use mechanical devices or other means to accelerate the defrosting process.
8. The freeze dryer contains flammable refrigerating equipment. Do not damage the refrigerant circuit. While handling, moving, and operating, avoid physical damage to reduce the risk of a leak.
9. Ensure that all accessories and pumps adhere to specifications.
10. Connect and disconnect the accessories only while the machine is unplugged.

! ATTENTION (may affect operational performance or service life)



1. The overall safe operation of the freeze dryer is the responsibility of the owner of the freeze dryer and their assigned operator(s), who in turn are responsible for ensuring the user manual guidance is applied, to ensure the safety and protection of personnel and the freeze dryer before, during, and after freeze drying operation.
2. The freeze dryer must be placed on a suitable surface of at least 77cm x 75cm (DxW), with a weight handling of minimum 100 kg. Allow an additional 20cm space in each dimension for ventilation. Allow good access to power connections in case of emergency. Ensure the freeze dryer can be unplugged easily during emergency.





3. Timely maintenance of the freeze dryer MUST be conducted to ensure continued safe operation and to optimize the freeze dryer's service life.
4. Only accredited and qualified professional repair technicians can open the freeze dryer or conduct required repairs. Persons performing repairs on the freeze dryer other than those selected or approved by the Company shall void any warranty for the product and are at risk of personal injury.
5. The freeze dryer must be used in a well-ventilated, dry environment.
6. The freeze dryer should be kept in a shaded area, without direct sunlight.
7. The freeze dryer must not be used in hazardous areas, or underground.

FIRE HAZARD!



1. R290 - Refrigerant used in this instrument.
2. H220 - Extremely flammable gas.
3. In the case of refrigerant leak, DO NOT operate the freeze dryer. Ensure the freeze dryer is disconnected from power supply and contact support immediately.

2.1 Panel Warning

 <p>Rear Panel</p>	<p>CAUTION – RISK OF FIRE OR EXPLOSION due to puncture of REFRIGERANT TUBING. Follow Handling Instructions Carefully. FLAMMABLE REFRIGERANT used.</p>
 <p>Door</p>	<p>WARNING – CRUSHING OF HANDS between the door and front panel, hinges, and door lock. Operate with caution.</p>
 <p>Chamber Wall, Shelves, and Trays</p>	<p>WARNING – FREEZING CONDITIONS. FROSTBITE HAZARD. Thermal gloves are required for extreme temperatures.</p>
 <p>Shelves and Trays</p>	<p>CAUTION – HOT SURFACE. Shelves and trays can be very hot. Thermal gloves are required for extreme temperatures.</p>

2.2 Lifting and Moving

To ensure safety when lifting or moving the freeze dryer, follow the below notes:

- Adjust the feet to reduce pinching hazard when positioning the freeze dryer.
- Lift or move the freeze dryer on the base only.
- Keep the freeze dryer oriented upward.
- Lift or move the freeze dryer after removing the shelf stack and disconnecting all accessories.
- Avoid dropping or any severe shock and vibration. Avoid sharp objects. Avoid penetration of the housing.
- Two people are required to lift. Lift from each corner of the freeze dryer, near the feet. Hold onto the base plate.

III. Sécurité

DANGER ! (peut causer de graves dommages aux biens ou faire des victimes)



1. Aucune pièce réparable par l'utilisateur ne se trouve à l'intérieur. L'entretien doit être effectué uniquement par du personnel autorisé.
2. Lisez ce manuel avant d'utiliser le lyophilisateur et vérifiez que toutes les connexions sont correctement fixées.
3. Exigences utilitaires : assurez-vous que les connexions électriques sont conformes aux normes locales et que l'alimentation électrique est compatible avec la plaque signalétique du lyophilisateur. La source d'alimentation doit être correctement mise à la terre.
4. N'utilisez pas le lyophilisateur si les câbles d'alimentation sont endommagés.
5. Le caisson et les étagères peuvent être extrêmement froids ou chauds. Vérifiez la température avant d'ouvrir le caisson ou de toucher ses composants internes, afin d'éviter les gelures ou les brûlures.
6. Ne forcez pas l'ouverture de la porte du caisson lorsque celle-ci est sous vide. Vous risqueriez de vous blesser gravement et d'endommager fortement le lyophilisateur.
7. N'exposez pas des parties du corps ou des corps entiers à la machine.
8. Le lyophilisateur a été exclusivement conçu pour la lyophilisation de produits solides ou liquides ininflammables, non corrosifs et non dangereux.
9. Il n'y a pas de batterie remplaçable par l'utilisateur.

AVERTISSEMENT ! (peut provoquer des dommages matériels ou corporels)



1. Lisez ce manuel avant d'utiliser le lyophilisateur.
2. Le lyophilisateur ne doit être ouvert que par le personnel d'entretien autorisé.
3. Ne placez pas d'objets lourds sur le lyophilisateur et ne vous tenez pas debout dessus.
4. N'installez aucun logiciel non autorisé sur le lyophilisateur.
5. N'approchez pas les mains de la porte du caisson pendant la descente de la pompe, afin d'éviter les blessures par écrasement.
6. Assurez-vous que toutes les ouvertures de ventilation ne sont pas obstruées et que l'air peut circuler librement.
7. N'utilisez pas de dispositifs mécaniques ou d'autres moyens pour accélérer le processus de dégivrage.
8. Le lyophilisateur contient des équipements frigorifiques inflammables. N'endommagez pas le circuit de réfrigération. Lors de la manipulation, du déplacement et de l'utilisation, évitez tout dommage physique afin de réduire le risque de fuite.
9. Veillez à ce que tous les accessoires et toutes les pompes soient conformes aux spécifications.
10. Branchez et débranchez les accessoires uniquement lorsque la machine est débranchée.

ATTENTION ! (conséquences possibles sur les performances opérationnelles ou la durée de vie)







1. La sécurité générale du fonctionnement du lyophilisateur relève de la responsabilité du propriétaire du lyophilisateur et de son (ses) opérateur(s) attitré(s), qui doivent à leur tour s'assurer que les conseils du manuel de l'utilisateur sont appliqués, afin de garantir la sécurité et la protection du personnel et du lyophilisateur avant, pendant et après l'opération de lyophilisation.
2. Le lyophilisateur doit être placé sur une surface appropriée d'au moins 77 cm x 75 cm (DxL), avec une capacité de charge d'au moins 100 kg. Prévoyez un espace supplémentaire de 20 cm pour la ventilation. Prévoyez un bon accès aux connexions électriques en cas d'urgence. Veillez à ce que le lyophilisateur puisse être débranché facilement en cas d'urgence.
3. La maintenance du lyophilisateur DOIT être effectuée en temps utile pour garantir un fonctionnement sûr et optimiser la durée de vie du lyophilisateur.
4. Seuls des techniciens de réparation professionnels accrédités et qualifiés peuvent ouvrir le lyophilisateur ou effectuer les réparations nécessaires. Les personnes effectuant des réparations sur le lyophilisateur autres que celles sélectionnées ou approuvées par la société annuleront toute garantie sur le produit et s'exposeront à des risques de dommages corporels.
5. Le lyophilisateur doit être utilisé dans un environnement sec et bien ventilé.
6. Le lyophilisateur doit être conservé dans un endroit ombragé, sans exposition directe au soleil.
7. Le lyophilisateur ne doit pas être utilisé dans les zones dangereuses ou souterraines.

RISQUE INCENDIE !



1. R290 - Réfrigérant utilisé dans cet appareil.
2. H220 - Gaz extrêmement inflammable.
3. En cas de fuite de réfrigérant, NE faites PAS fonctionner le lyophilisateur. Assurez-vous que le lyophilisateur est débranché de l'alimentation électrique et contactez immédiatement le service d'assistance.

3.1 Avertissements du panneau

 <p>Panneau arrière</p>	<p>ATTENTION -RISQUE D'INCENDIE OU D'EXPLOSION dû à la perforation du TUBE RÉFRIGÉRANT. Suivez attentivement les instructions de manipulation. RÉFRIGÉRANT INFLAMMABLE utilisé.</p>
 <p>Porte</p>	<p>AVERTISSEMENT -ÉCRASEMENT DES MAINS entre la porte et le panneau avant, les charnières et la serrure de la porte. Utilisez avec précaution.</p>
 <p>Mur, étagères et plateaux du caisson</p>	<p>AVERTISSEMENT -CONDITIONS DE GEL. RISQUE DE GELURES. Des gants thermiques sont nécessaires en cas de températures extrêmes.</p>
 <p>Étagères et plateaux</p>	<p>ATTENTION -SURFACE CHAUDE. Les étagères et les plateaux peuvent être très chauds. Des gants thermiques sont nécessaires pour les températures extrêmes</p>

3.2 Levage et déplacement

Pour garantir la sécurité lors du levage ou du déplacement du lyophilisateur, veuillez respecter les consignes suivantes:

- Ajustez les pieds pour réduire le risque de pincement lors du positionnement du lyophilisateur.
- Soulevez ou déplacez le lyophilisateur sur la base uniquement.
- Maintenez le lyophilisateur orienté vers le haut.
- Soulevez ou déplacez le lyophilisateur après avoir retiré la pile d'étagères et débranché tous les accessoires.
- Évitez les chutes, les chocs et les vibrations. Évitez les objets pointus. Évitez toute pénétration dans le boîtier.
- Deux personnes sont nécessaires pour soulever l'appareil. Soulevez chaque coin du lyophilisateur, près des pieds. Tenez la plaque de base.

IV. Introduction

Users are advised to carefully read this manual prior to operating the freeze dryer so that they are aware of all precautions outlined and to ensure operation is in accordance with the instructions contained within this manual.

4.1 After Sales Support

If problems are encountered or technical support is required when installing or using the instrument, please contact sales@hollandgreenscience.eu

The company may provide technical assistance and information regarding the instrument or equipment or service without charge at its sole discretion. Buyer assumes sole responsibility for any reliance on or use of such assistance and information, and the company makes no warranty thereon.

Upon contact the following information is required:

- Product serial number (located on the instrument nameplate)
- Description of issue or problem
- Method and or operating steps you have undertaken towards resolution.
- Your contact details inclusive of telephone number and email address.

4.2 Proper Use

The freeze dryer has been exclusively designed for freeze drying of non-flammable, non-corrosive and non-hazardous solids or liquid containing products. Freeze drying flammable, corrosive or hazard solids or liquids could damage the freeze dryer, void the warranty and cause accidents. The freeze dryer must not be used in an environment with a potentially hazardous or flammable atmosphere. The freeze dryer must not be used underground.

The freeze dryer is primarily designed for non-residential commercial use and is to be used only in conjunction with accessories recommended within this manual.

4.3 Product Disposal

The freeze dryer should not be disposed of with other household or electronics waste. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle responsibly to promote the sustainable reuse of material resources. To return your used device, use the return and collection systems, or contact the distributor where the product was purchased. This product contains polyurethane (PU) and R290 (propane). PU is flammable and recyclable. R290 is an extremely flammable refrigerant.

4.4 Warranty

The warranty period for the Xiros Mikro freeze dryer is one (1) year from the date of

shipping of the product. Please note that when a unit is replaced, the warranty does not restart. Instead, the warranty for the replacement unit transfers from the original, covering only the remaining balance of the original warranty period. Please ensure that you are familiar with the warranty period for the unit.

4.5 Customer Responsibilities & Warranty Conditions

To maintain your warranty coverage, the consumer is responsible for the following:

- **Inspection Upon Delivery:** As soon as your shipment arrives, please inspect the freeze dryer for any damage that could have occurred during transit. Please take care when removing the unit from the crate and save the original packaging until you have turned it on and registered your Xiros Mikro small freeze dryer.
- **Proper Use:** Operating the appliance in accordance with the official instructions provided in the user manual included with your Xiros Mikro freeze dryer.
- **Proper Installation:** Installing the unit according to the instruction manual and in compliance with all local electrical codes.
- **Electrical Supply:** Connecting the unit to a dedicated and properly grounded power supply of sufficient voltage. This includes responsibility for replacement of blown fuses, repair of loose connections, or addressing any defects in home wiring.
- **Operating Environment:** Placing the unit in a clean, open area with ample airflow, where the ambient temperature remains between 10°C and 30°C.
- **Handling & Transport:** Coordinating delivery and moving the appliance from the delivery truck to its final location. **Note:** Freight drivers typically deliver to the curb side or end of the driveway; it is the customer's responsibility to arrange for assistance moving the unit inside if needed. To prevent damage, never lift the unit by the door.
- **Routine Maintenance:** Properly maintaining the freeze dryer and vacuum pump according to the detailed schedule in the owner's manual.
- **Consumable Parts:** Performing minor part replacements, maintenance, and repairs as directed by Xiros Mikro support. This includes, but is not limited to, O-rings, gaskets, seals, valves, and filters.
- **Parts Return:** Returning any requested failed parts to Holland Green Science for evaluation when a prepaid shipping label is provided.
- **Delivery Access:** Ensuring the unit is delivered to a location where the freight truck can safely and legally access. It is your responsibility to communicate any access restrictions (e.g., narrow roads, low-hanging wires, steep inclines) to the freight company in advance.



V. Technical Specifications

Product Name	Xiros Mikro
Product SKU	20107002C
Product Weight	87 kg
Max Vacuum Pump Current	AC120V, 50/60Hz, 2.8A AC220-240V, 50/60Hz, 1.6A (Not for UL-mark)
Total Condenser Volume	11 liters
Ice Condenser Capacity	8 kg
Condenser Performance	4 kg per 24 hours
Ice Condenser Temperature	< -45 °C
Ultimate Vacuum	1.5 x 10 ⁻¹ mbar
Minimum Shelf Temperature	-40 °C
Maximum Shelf Temperature	+60 °C
Tray Dimensions (mm)	9-shelf tray: 200 (W) x 450 (L) x 12 (H) 7-shelf tray: 200 (W) x 450 (L) x 15 (H) 5-shelf tray: 200 (W) x 450 (L) x 20 (H) 3-shelf tray: 200 (W) x 450 (L) x 20 (H)
Number of Trays	5
Shelf Stack Options	9, 7, 5, 3 shelf
Shelf Stack Dimensions (mm)	465 (L) x 210 (W) x 240 (H)
Typical Product Capacity	9.6 kg at 80% moisture
Distance Between Trays	5-shelf: 33 mm
External Dimensions (mm)	770 (D) x 748 (W) x 507 (H)
Operating Altitude	< 2000m
Environmental Specification	10°C to 30°C 10% to 80%, non-condensing
Area of Operation	Indoors only
Overvoltage Category	Category II
Voltages and Frequencies	For North America with UL mark AC120V 50/60Hz 12A For the rest of the world AC100-120V, 50/60Hz 12A AC220-240V, 50/60Hz 8A
Pollution Degree	2

VI. Precautions

- The freeze dryer should be connected to electrical outlets safely. Operate with dry hands to minimize the risk of electric shock.
- A suitably fused power supply (according to the freeze dryer's nameplate) is required. An alternative supply voltage or frequency may damage the freeze dryer or reduce its working lifespan.
- The electrical power cord is designed to be connected to the power outlet without knots, sharp bends, or heavy materials placed on the cord.
- If the power cable is damaged, replace with same specification or contact sales@hollandgreenscience.eu.
- The freeze dryer must be properly grounded according to local electrical codes.
- No residual water or foreign matter should be present inside the chamber prior to freeze drying.
- If an abnormal sound, excessive heat, smoke, etc. is detected, stop the process immediately, disconnect the freeze dryer from the power supply, and contact sales@hollandgreenscience.eu. The freeze dryer operates with a low level of noise when running; if any significant changes occur, stop and contact support.
- If a power outage occurs when operating the freeze dryer, open the drain valve and let the chamber pressure return to normal before opening the door to retrieve the product.
- When the freeze-drying process has finished, first turn off the vacuum pump, and then open the drain valve.
- Ensure the door seal ring and the chamber door remain clean. Only clean with soapy water and avoid using solvents or other cleaning agents.
- Do not rapidly power on/off the freeze dryer. Wait for at least a minute after powering off the freeze dryer to power it on again.
- It is recommended that the freeze dryer be unplugged from the wall outlet, or the outlet be switched off, when not in use.
- Water may escape from the freeze dryer if it is operated outside of recommended operating conditions. Regularly check the vacuum pump for potential water contamination. Depending on the pump, replace the oil, and service as needed.
- Do not rapidly enable the compressor, as it may damage the component. Wait at least 30 seconds before enabling it again.
- While the compressor is running, do not directly power off the freeze dryer. Turn it off from the software first, before turning off the power.
- When not in use, keep the drain valve and the chamber door open to prevent micro bacterial and mold growth due to residual moisture. Only seal off for storage if the chamber is entirely dry.

VII. Components Overview





Power Socket

Fuse (250VAC F15AL)

Vacuum Pump Supply Output



Drain Valve Closed



Drain Valve Open

VIII. Software Operation

Immediately after powering on the freeze dryer, the home screen is shown (Figure 1).




Figure 1

To go back to the home screen, tap the “Home” icon that is present on most screens (Figure 2).




Figure 2

8.1 AutoDry

AutoDry automatically dries the product using an intelligent dryness detection algorithm. AutoDry cannot guarantee complete product dryness for each run performed. If the dried product is found to be too moist, proceed to “8.1.2 Editing AutoDry Parameters”.

From the Home screen, tap on “AutoDry” (Figure 3).




Figure 3

Figure 4 shows the AutoDry Start screen. Tap on “Start” to begin the freeze-drying process.




Figure 4

If in operation (i.e. Program or Defrost is running), Figure 5 shows an error message with the “Start” button blocked.




Figure 5

To proceed, stop Defrosting or the Program by navigating onto the respective screen, and tap on “Stop” or “Cancel”.

8.1.1 AutoDry Process

AutoDry is split into 4 stages. Chronologically, they are Freezing, Evacuation, Drying, and Storage. A typical run will last for roughly 14-30 hours, depending on the moisture content of the product, and the parameters specified for the AutoDry algorithm.

Figure 6 shows the first AutoDry stage: the Freezing stage.




Figure 6

"Estimated Running Time" is displayed at the bottom of the screen. This time is based on the average runtime of previous successful AutoDry runs. Therefore, the freeze dryer first needs to almost complete a run before the displayed time is self-corrected. Adjusting any of the AutoDry parameters will reset this timer.

If there has been a power outage during a run, timing information will not be recorded for that run, and the remaining time may become inaccurate.

This timer works best using consistent product type, size and weight.

Figure 7 shows the Evacuation stage.




Figure 7

Figure 8 shows the Drying stage.




Figure 8

Figure 9 shows the Storage stage.




Figure 9

8.1.2 Editing AutoDry Parameters

To begin editing the parameters, tap on “Edit Parameters” from the AutoDry Start screen (Figure 10).




Figure 10

Figure 11 shows the Parameter Editor screen. This screen provides significant control over the behavior of the AutoDry algorithm. Only 3 of the 4 stages can be modified with AutoDry: Freezing, Drying, and Storage.




Figure 11

“Freeze Timer Start” (Figure 12) specifies the temperature the shelves must reach before the timer begins. “Freeze Timer Duration” specifies how long the shelves must be below the “Freeze Timer Start” temperature. This helps to ensure that the product is thoroughly frozen before drying starts.




Figure 12

Figure 13 shows the Drying stage. You may increase the maximum drying temperature in conjunction with reducing sensitivity and adding more drying time. Higher temperatures and lower sensitivity will achieve shorter drying times but may have negative impacts on product quality. Reducing drying temperatures and increasing sensitivity will likely produce better quality results.




Figure 13

The Storage stage (Figure 14) specifies the freeze dryer operation at the conclusion of a run. We recommend storage above ambient temperature, to prevent rehydration.




Figure 14

To restore to the default parameters for AutoDry, tap on “Reset”. The following popup will display (Figure 15):




Figure 15

Tapping on “Confirm” will overwrite all previous changes to the AutoDry Parameters screen.

Note: Adjusting any parameters or resetting the values will require AutoDry to re-learn the duration of the run, and therefore will not be able to provide the duration hint until it has nearly completed a run.

8.2 Program

The freeze dryer has up to 12 customizable programs available, each with 16 steps. From the Home screen, tap on “Program” (Figure 16) to access a list of available programs (Figure).




Figure 16

To review the remaining 6 available programs, tap on the right arrow (Figure 17).




Figure 17

Figure 18 shows the remaining 6 programs. To go back, tap on the left arrow.



Figure 18

To edit one of these programs, tap on one of the programs (Figure 19) to enter the program editor. This will be explored in the later chapters.



Figure 19

8.2.1 AutoDry/Program Stages

The freeze dryer Programs and AutoDry are broken down into 4 stages. Chronologically, they are as follows: Freezing, Evacuation, Drying, and Storage. See section “**8.2.2 Writing a Program & Parameters Explained**” for more details on each of the parameters, and how they control each stage.

- Freezing: This stage is used to cool the shelves, trays and product. This stage generally lasts for a few hours and will extend with increased load.
- Evacuation: This stage evacuates the pressure inside the chamber before applying energy. This stage generally lasts for a few minutes.
- Drying: This stage involves the sublimation process where energy (heat) is transferred to the product under low chamber pressure.
- Storage: This is the final stage of the freeze-drying process. This determines a safe pressure and temperature to hold the product at until it gets removed from the freeze dryer.

8.2.2 Writing a Program & Parameters Explained

From the Program Selector screen discussed in chapter 8.2, tap on “Program 1” to enter the Program Editor screen, as shown in Figure 20.




Figure 20

In Figure 20 the entries are only an example, as the parameters are dependent on the product. VO is explored in “8.2.6 Program Vacuum Oven (VO)”.

8.2.2.1 Freezing Stage

The Freezing stage (Figure 21) is controlled by the “Freeze Timer Start” and “Freeze Timer Duration”. “Freeze Timer Start” specifies at what shelf temperature the timer will start, and “Freeze Timer Duration” specifies how long the shelves must stay under the specified temperature before proceeding to the Evacuation stage.





Figure 21

8.2.2.2 Evacuation Stage


The Evacuation stage (Figure 22) is controlled by Step 1's pressure. Once the specified pressure is reached, the freeze dryer will proceed to the next stage.



Freeze Timer Start (°C)		Freeze Timer Duration (min)		
-20.0		240		
Drying				
Drying Step	Temperature (°C)	Time (min)	Pressure (mbar)	Incl. Ramp Time
Step 1	-30.0	120	1.00	<input type="checkbox"/>
Step 2	-20.0	120	0.80	<input type="checkbox"/>
Step 3	-10.0	120	0.60	<input type="checkbox"/>
Step 4	0.0	120	0.40	<input type="checkbox"/>
Step 5	10.0	120	0.30	<input type="checkbox"/>
Step 6	20.0	120	0.30	<input type="checkbox"/>
Step 7	30.0	120	0.30	<input type="checkbox"/>
Step 8	0.0	0	0.00	<input type="checkbox"/>

Figure 22

The Evacuation stage will only run for up to 15 minutes. If the run has exceeded that time and the pressure is more than 0.5 mbar above from the setpoint, the freeze dryer will abort the run and enable the compressor in the manual page. The below message will appear, indicating that the Evacuation stage has failed (Figure 23).



EVACUATION FAILED

Xiros Mikro cannot evacuate.
The operation is aborted.

Pressure reached (mbar)	10.53
Pressure goal (mbar)	1.50

- ✓ Hose clamps are secure
- ✓ Drain valve is closed
- ✓ Chamber has recently defrosted
- ✓ Vacuum pump is switched on

The compressor has remained on to protect the product

Confirm

Figure 23

If VO is enabled (see “8.2.6 Program Vacuum Oven (VO)”), then the compressor will not be enabled. Return to the Home screen by tapping on “Confirm”.

8.2.2.3 Drying Stage

In the Drying stage (Figure 24), up to 16 drying steps can be specified. They run sequentially, and are composed of a temperature and pressure setpoint, time, and the option of “Include Ramp Time”.

Incl. Ramp Time: This controls when the timer will start for the Drying stage. When the checkbox is checked, the timer will immediately start for the Drying stage. When all steps have this option checked, the total process time equals the sum of all drying steps (independent from how long it takes the freeze dryer to ‘ramp’ from step to step). When the checkbox is unchecked, the timer will only start once the shelves reach the temperature setpoint.

Freeze Timer Start (°C)		Freeze Timer Duration (min)		
-20.0		240		
Drying				
Drying Step	Temperature (°C)	Time (min)	Pressure (mbar)	Incl. Ramp Time
Step 1	-30.0	120	1.00	<input type="checkbox"/>
Step 2	-20.0	120	0.80	<input type="checkbox"/>
Step 3	-10.0	120	0.60	<input type="checkbox"/>
Step 4	0.0	120	0.40	<input type="checkbox"/>
Step 5	10.0	120	0.30	<input type="checkbox"/>
Step 6	20.0	120	0.30	<input type="checkbox"/>
Step 7	30.0	120	0.30	<input type="checkbox"/>
Step 8	0.0	0	0.00	<input type="checkbox"/>

Figure 24

Only the first 8 steps are visible on the first page. Tap on the right arrow to navigate to the next screen, which contains the remaining 8 steps (Figure 25).

Storage Temperature (°C)		Storage Pressure (mbar)		
30.0		3.00		
Drying				
Drying Step	Temperature (°C)	Time (min)	Pressure (mbar)	Incl. Ramp Time
Step 9	0.0	0	0.00	<input type="checkbox"/>
Step 10	0.0	0	0.00	<input type="checkbox"/>
Step 11	0.0	0	0.00	<input type="checkbox"/>
Step 12	0.0	0	0.00	<input type="checkbox"/>
Step 13	0.0	0	0.00	<input type="checkbox"/>
Step 14	0.0	0	0.00	<input type="checkbox"/>
Step 15	0.0	0	0.00	<input type="checkbox"/>
Step 16	0.0	0	0.00	<input type="checkbox"/>

Figure 25

8.2.2.4 Storage stage

The freeze dryer will maintain a specified temperature and pressure that the product will be left at once the run has finished (Figure 26). Specify a pressure such as 1 mbar to avoid additional drying whilst in Storage and specify a temperature between 25-40C to prevent the product from absorbing moisture when it is reintroduced to the atmosphere.

The screenshot shows the 'Xiros Mikro' interface for 'Program-1'. It features a 'Drying' table with columns for 'Drying Step', 'Temperature (°C)', 'Time (min)', 'Pressure (mbar)', and 'Incl. Ramp Time'. Below the table, a pink highlighted box contains the 'Storage Temperature (°C)' set to 30.0 and 'Storage Pressure (mbar)' set to 3.00. Navigation arrows are visible on the left and right sides of the interface.

Drying				
Drying Step	Temperature (°C)	Time (min)	Pressure (mbar)	Incl. Ramp Time
Step 9	0.0	0	0.00	<input type="checkbox"/>
Step 10	0.0	0	0.00	<input type="checkbox"/>
Step 11	0.0	0	0.00	<input type="checkbox"/>
Step 12	0.0	0	0.00	<input type="checkbox"/>
Step 13	0.0	0	0.00	<input type="checkbox"/>
Step 14	0.0	0	0.00	<input type="checkbox"/>
Step 15	0.0	0	0.00	<input type="checkbox"/>
Step 16	0.0	0	0.00	<input type="checkbox"/>

Storage Temperature (°C)	Storage Pressure (mbar)
30.0	3.00

Figure 26

8.2.3 Renaming a Program

The freeze dryer provides functionality to rename programs. On either the first or second Program Editor screen, tap on the program's name to prompt the keyboard function (Figure 27).




Figure 27

Type in a name, and then tap “Enter” to confirm the input. The name will update. The program name will be displayed on the Program Selection, Program Progress, and System Idle screens.

8.2.4 Loading a Preset Program

The freeze dryer comes with 4 preset programs as examples. To load a program, first select a program to load a preset program from the Program Selector screen. The button to load programs is located on the top-right of the Program Editor (Figure 28).




Figure 28

To prompt the UI for selecting a loadable program, tap on “Load” (Figure 29).




Figure 29

The following popup displays all the available programs to load (Figure 30). SafeDry, Liquid 2L, Candy, and Botanical example programs are available. Tapping on one of them will show a popup with product preparation instructions.




Figure 30

Tap on "SafeDry". Figure 31 is the popup that appears. Tapping on "Load" will close the popup and write the SafeDry program onto the selected program.




Figure 31


8.2.5 Starting a Program

Before starting the program, ensure the following:

- Product is evenly distributed between each tray for optimal drying.
- The drain valve is closed.
- The chamber is defrosted.
- The chamber door contact surface is clear of ice, and the door is closed.

It is important to defrost the freeze dryer before conducting a run. In addition, occasionally opening the ballast valve on the pump to release moisture buildup is recommended.

To run a program, while in the program editor, tap on the right arrow (Figure 32) to navigate to the Program Start screen (Figure 33).






Figure 33

When you navigate to the Program Start page, it may take a few seconds to initialize. Wait until it initializes before attempting to start. If under an operation (e.g. Defrosting or AutoDry), you will get the following screen (Figure 34).




Figure 34

You will need to stop the other mode (such as Defrost) to start a program. If no other mode is running, and it has successfully initialized, you may tap on "Start" to begin the program (Figure 35).




Figure 35

Verify the program is running by the sound of the compressor starting and the Program Progress screen appearing (Figure 36).




Figure 36

The Program Progress screen displays real-time values of the operation of the current program (Figure 36). Some of the readings are as follows:

Elapsed Time: The number of minutes that have passed since the program started.

Remaining Section Time: The timer for how much time is left on the current stage/drying step.

Remaining Total Time: The cumulative total of all the times from each stage and drying steps. This is not the estimated remaining time, but rather the remaining timer time.

Note: Most timers don't start unless the preconditions are met. For example, if the "Incl. Ramp Time" is unchecked, the drying step timer starts once the shelf temperatures are reached. This is indicated with a message under the "Cancel" button.

The Program Progress screen also displays the current stage. There are four stages in a Program "Freezing" (Figure 36), "Evacuating" (Figure 37), "Drying" (Figure 38) and "Storage" (Figure 39).




Figure 37

Note: If the pressure isn't reached within a certain amount of time, the evacuation stage is aborted. This is explored in "8.2.2.2 Evacuation Stage".




Figure 38




Figure 39

When the Stage/Drying step displays “Storage” (Figure 39), the Program is finished and will maintain temperatures and pressure based on the storage parameters. Stopping a program is explored in “8.3 Extra Drying”.

8.2.6 Program Vacuum Oven (VO)

Vacuum Oven (VO) provides the ability for the freeze dryer to serve as a vacuum oven. Enabling VO skips the Freezing stage, and the compressor remains off throughout the run.

VO is enabled on the Program Editor’s first screen, by checking the VO box on the top left (Figure 40).




Figure 40

Once VO Mode is enabled, the Freezing stage is no longer editable and is skipped (Figure 41).




Figure 41

8.3 Extra Drying

Once a program or AutoDry run is in storage, the operator will have the ability to verify that their product is dry before ending the program. After tapping “Stop” on the Program Progress screen, a popup will appear (Figure 42).




Figure 42

While this popup is shown, the vacuum control valve will remain closed. This allows the operator to open the drain valve, retrieve the product and test it, while the shelves and compressor remains on. Tapping on “Yes” will stop the run, while tapping on “Cancel” will close the popup and re-enable the opening of the vacuum control valve.

To start the extra drying, tap on “No”, and the following popup will appear (Figure 43):




Figure 43

Figure 43 shows the Extra Drying Parameters popup. This popup will allow you to specify an extra drying step for your Program. Specifying pressure is not available in AutoDry (Figure 44).




Figure 44

The parameters for the extra drying will default, depending on whether the extra drying is being applied to AutoDry or Program Mode. If AutoDry, the Temp and Time parameters will be based on “Max Drying Temp” and “More Drying Time” respectively, while for Program, the parameters will be based on the last step for the Program.

After tapping on “Confirm”, the popup will disappear, and “Extra Drying” will be visible on the Drying Stage section of the Program Progress screen (Figure 45).




Figure 45

8.4 Power Recovery

If a transient power outage occurs during a Program or AutoDry run, the freeze dryer will automatically recover back to the last state before the freeze dryer lost power. A popup will appear as soon as power is restored to the freeze dryer (Figure 46).




Figure 46

This is designed to be completely autonomous. The freeze dryer will continue the drying process irrespective of the duration or frequency of the power outages. It is up to the operator to determine the health of the product.

During an AutoDry run, recovering from power loss will prevent timing information from being recorded for that run, and the remaining time may become inaccurate.

8.5 Manual Control

To enter the Manual Control screen, tap on the “Manual” button (Figure 47).




Figure 47

Figure 48 shows the Manual Control screen.




Figure 48

On the right-hand side (Figure 49), temperature and pressure setpoints can be specified for the freeze dryer to maintain (if enabled).




Figure 49

The left-hand side provides relay control of the freeze dryer. They are as follows:

- Compressor
- Vacuum Valve
- Vacuum Pump
- Shelf Heaters (upper & lower)

Tap on those relays (Figure 50) to enable them.




Figure 50

Note: For the freeze dryer to maintain the pressure setpoint, enable both the Vacuum Pump and Vacuum Valve. The Vacuum Valve will not open if the pump isn't enabled, and a warning will display (Figure 51). This is to prevent suck-back contamination from the pump.




Figure 51

Ensure that there are no liquids, and your product is fully frozen, before evacuating. Evacuating while there is still liquid in the chamber may damage or cause premature wear to the vacuum pump.

8.5.1 Defrost Mode

While the freeze dryer is in an idle state, tap on “Defrost” from the home screen to enter the Defrost screen (Figure 52).




Figure 52

Figure 53: The Defrost screen prior to starting.




Figure 53

Before starting the Defrost process, confirm:

- Product has been taken out.
- Drain valve is open.
- Chamber door is closed.

Set the Defrost Temperature and Defrost Time as desired (Figure 54). A 50°C defrost temperature and 2 hours hold time are generally adequate.




Figure 54

Tap on “Start” (Figure 55).




Figure 55

The “Defrosting” icon will appear (Figure 56).




Figure 56

The System Defrost screen will display the real-time values for the lower and upper shelf temperature and chamber temperature.

Tap “Cancel” to stop the defrosting process (Figure 57).




Figure 57

8.6 System Settings

Tap "Settings" on the home screen (Figure 58).




Figure 58

(Figure 59). Rebooting the display is not the same as power cycling the freeze dryer. After updating, perform a power cycle using the main power switch located under the display.




Figure 59

The screen timeout specifies the time it takes to display the System State since the last touch interaction. Setting it to 0 minutes will disable it. The System State screen (Figure 60) shows a general overview of the current state of the freeze dryer.




Figure 60

8.7 History Screen

The History screen shows data from the freeze dryer over time, in 15 second intervals. To access this screen, tap on “History” from the home screen (Figure 61).




Figure 61


Figure 62 is the History screen. Tap on “Page Up” and “Page Down” to navigate through the history of the temperature/pressure sensors with the associated time stamps.

Time	Date	Upper Shelf Temp (°C)	Lower Shelf Temp (°C)	Condenser Temp (°C)	Pressure (mbar)
17:02:51	2023/06/21	16.6	16.3	16.5	99.99
17:03:06	2023/06/21	16.7	16.3	16.5	99.99
17:03:21	2023/06/21	16.6	16.2	16.5	99.99
17:03:36	2023/06/21	16.6	16.2	16.5	99.99
17:03:51	2023/06/21	16.5	16.3	16.6	99.99
17:04:06	2023/06/21	16.6	16.3	16.6	99.99
17:04:21	2023/06/21	16.6	16.3	16.5	99.99
17:04:36	2023/06/21	16.5	16.3	16.5	99.99
17:04:51	2023/06/21	16.5	16.3	16.5	99.99
17:05:06	2023/06/22	16.4	16.2	16.5	99.99

Figure 62

8.7.1 History Export

To export the History of the freeze dryer, insert a USB formatted as FAT32 into the USB-A port of the freeze dryer. After inserting the USB, navigate to the History screen and tap on “Export” (Figure 63).



The screenshot shows the 'Xiros Mikro' interface with a 'HISTORY' screen. At the top, there are three navigation buttons: 'Page Up', 'Page Down', and 'Export'. A hand icon is shown tapping the 'Export' button. Below the buttons is a table with the following data:

Time	Date	Upper Shelf Temp (°C)	Lower Shelf Temp (°C)	Condenser Temp (°C)	
17:02:51	2023/06/21	16.6	16.3	16.5	99.99
17:03:06	2023/06/21	16.7	16.3	16.5	99.99
17:03:21	2023/06/21	16.6	16.2	16.5	99.99
17:03:36	2023/06/21	16.6	16.2	16.5	99.99
17:03:51	2023/06/21	16.5	16.3	16.6	99.99
17:04:06	2023/06/21	16.6	16.3	16.6	99.99
17:04:21	2023/06/21	16.6	16.3	16.5	99.99
17:04:36	2023/06/21	16.5	16.3	16.5	99.99
17:04:51	2023/06/21	16.5	16.3	16.5	99.99
17:05:06	2023/06/22	16.4	16.2	16.5	99.99

Figure 63

After tapping on “Export”, a dialogue box will appear showing “Exporting Data...”. Wait until “Exporting Complete” is shown and then disappears, then remove the USB. Insert the USB into a computer. The exported file will be called “History-YYYY-MM-DD-HH-mm-SS.xls”, located onto the root directory of the USB. The filename specifies the timestamp of when the history was exported.

IX. Fault Diagnosis

Before conducting a diagnosis, verify that the chamber is clean, dry, and empty of product.

Error	Solutions
Ultimate Vacuum of 1.5×10^{-1} mbar not reached.	Check that the vacuum pump is correctly connected to the main unit with the clamp (tightly).
	Check that the drain valve is closed and is installed hand tight.
	Check that the vacuum hose and rings are installed correctly. Verify that there is no damage on the rubber rings, and that it is clear from particles.
	Check that the vacuum pump functions correctly, and that the vacuum pump oil is clear (if oil sealed vacuum pumps are used).
	Verify that the ballast valve is in a closed position on the pump (if fitted).
	Verify that the freeze dryer is defrosted, and that the inner door ring contacts a clean surface.
	Check that the compressor is enabled and has been running for at least 1 hour.
	If the error remains, contact support.
The chamber temperature is not cooling to specification.	Let the compressor run for at least 1 hour.
	Verify that no significant sounds are coming from the compressor. The freeze dryer should operate quietly. Contact support if the freeze dryer makes significant or unfamiliar noises.
	Check that shelves are not enabled in the Manual screen.
	Make sure the compressor is running by listening for sound.
	If compressor is "ON" but the temperature inside the chamber is still high, check the housing vents for blockage and dust. Remove blocking items and clean as required. If substantial dust is visible on the condenser inside the housing, contact support.
	Ensure that there is enough clearance around the vents, and that ambient temperature of the room is below 25°C.
The unit is making a strange sound.	Turn off the machine and contact support immediately.
The clock resets after each power cycle.	Contact support.

X. Cleaning and Maintenance

10.1 Cleaning

- Water from the product will collect inside of the chamber during drying. Be sure to defrost after each run and clean the inside of the chamber periodically by wiping down the chamber walls with washing detergent.
- Clean the trays after every freeze-drying process with mild, non-corrosive detergent. Wipe dry afterwards.
- Periodically check the chamber door seal ring and contact surface is clean from contaminants, to ensure performant evacuation.

10.2 Maintenance

- Check the chamber door seal ring regularly to make sure it is not damaged or worn.
- If an oil vacuum pump is used, check the vacuum oil after every 100 hours of operation, and replace the vacuum pump oil if the level is low or if the oil is dirty.
- If a dry pump is used, refer to its User's Manual.
- Clean instrument vents regularly to prevent dust build up.