

Holland
Green
Science



20209001

Aeolus

Vacuum Controller

User Manual

Table of Contents

I. LIMITED WARRANTY & DISCLAIMER.....	3
II. Safety	6
III. Introduction.....	7
IV. Technical Specifications	8
V. Inspection.....	9
VI. Operating Instructions	10
VII. Fault Diagnosis.	16
VIII. Maintenance and Cleaning.....	17

I. LIMITED WARRANTY & DISCLAIMER

- The intellectual property rights for this manual belongs to Holland Green Science Corporation (hereinafter referred to as “the Company”).
- Products of the company are under the patent protection of China and other countries and regions (inclusive of obtained patents and patents currently under application).
- The company reserves the right to change product specifications without notice.
- The contents of the manual are subject to change without notice.
- The contents of the manual must not be copied, extracted, or modified in any form without the Company’s permission.

LIMITED WARRANTY & DISCLAIMER

Products sold by the Company are intended solely for the use(s) specified by the manufacturer and are not to be used for any other purposes including, but not limited to, unauthorized commercial purposes. Buyer warrants and represents that Buyer will properly test, use, qualify, and/or validate each Product for intended use, and, to the extent authorized, manufacture and market any final articles made from the Company’s Products in accordance with the practices of a reasonable person who is an expert in the field, and in strict compliance with all applicable national, state and local laws and regulations. Buyer acknowledges and understands that the purchase of Products only conveys to Buyer the non-transferable right for only the Buyer to use the Products purchased in compliance with the applicable manufacturer’s Published User Statement, Limited User Statement, or Limited License, if any. While this manual contains certain representations of the potential dangers from the instrument or equipment, Buyer is solely responsible for conducting any research necessary to learn the hazards involved for any of its intended uses of the instrument or equipment purchased from the Company and to properly warn its customers, employees and other personnel who may be exposed to such instrument or equipment of any risks involved in using or handling the instrument or equipment.

The Company warrants that this product will operate or perform substantially in conformance with the Company’s published specifications and be free from defects in material and workmanship, when subjected to normal, proper and intended storage and usage by properly trained personnel, for the duration of the manufacturer’s warranty period set forth in this product documentation, published specifications or package inserts. If a manufacturer’s warranty period is not specified in the product documentation, published specifications or package inserts, the warranty period shall be one (1) year from the date of Bill of Sale to Buyer for equipment, and the earlier of any stated expiration dates or ninety (90) days for all other products (the “Warranty Period”). The Company agrees during the Warranty Period, to repair or replace, at the Company’s option, defective Products so as to cause the same to operate in substantial conformance with said published specifications; provided that Buyer shall (a) promptly notify the Company in writing upon the discovery of any defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) after the Company’s review, the Company will provide Buyer with service data and/or a Return Material Authorization (“RMA”), which may include biohazard decontamination procedures and other product-specific handling instructions, then, if applicable, Buyer may return the defective Products to the Company

with all costs prepaid by Buyer. Replacement parts may be new or refurbished, at the election of the Company. All replaced parts shall become the property of the Company. Shipment to Buyer of repaired or replacement Products shall be made in accordance with the Delivery provisions of the Terms & Conditions agreed to by the purchase of this product. Consumables are expressly excluded from this warranty.

Notwithstanding the foregoing, Products supplied by the Company but not manufactured by the Company are not warranted by the Company, and the Company disclaims all warranties associated with the accuracy of any original manufacturer or third-party supplier Product information, but the Company agrees to assign to Buyer any warranty rights in such Product that the Company may have from the original manufacturer or third-party supplier, to the extent such assignment is allowed by such original manufacturer or third-party supplier.

In no event shall the Company have any obligation to make repairs, replacements or corrections required, in whole or in part, as the result of (i) normal wear and tear, (ii) accident, disaster or event of force majeure, (iii) misuse, fault or negligence of or by Buyer, (iv) use of the Product in a manner for which they were not designed, (v) causes external to the Product such as, but not limited to, power failure or electrical power surges, (vi) improper storage and handling of the Product or (vii) use of the Product in combination with equipment or software not supplied by the Company. If the Company determines that Product for which Buyer has requested warranty services are not covered by the warranty hereunder, Buyer shall pay or reimburse the Company for all costs of investigating and responding to such request at the Company's then prevailing time and materials rates. If the Company provides repair services or replacement parts that are not covered by this warranty, Buyer shall pay the Company therefor at the Company's then prevailing time and materials rates.

Any installation, maintenance, repair, service, relocation or alteration to or of, or other tampering with, the product performed by any person or entity other than the company without the company's prior written approval, or any use of replacement parts not supplied by the company, shall immediately void and cancel all warranties with respect to the product.

The obligations created by this warranty statement to repair or replace a defective product shall be the sole remedy of buyer in the event of a defective product. Except as expressly provided in this warranty statement, the company disclaims all other warranties, whether express or implied, oral or written, with respect to the products, including without limitation all implied warranties of merchantability or fitness for any particular purpose or use. The company does not warrant that the products are error-free or will accomplish any particular result.

Nothing in this manual or in the terms and agreement agreed to at the time of purchase shall be construed to provide for the transfer, assignment or alienation of the limited warranties provided by the company herein to any other party other than buyer, including but not limited to, successors, heirs, subsidiaries or any other third person.

The warranties identified above are the sole and exclusive warranties with respect to this product and are in lieu of all other warranties, express or implied. All other warranties are expressly disclaimed, including, without limitation, those of non-infringement, merchantability, or fitness for a specific purpose or use (whether arising from statute, or otherwise in law, or from a course of dealing, or usage of trade). This warranty is not transferable from the original purchaser to a subsequent owner. further, the company is not liable in cases of deliberate, negligent or accidental misuse of the products, use with inappropriate reagents or consumables, damage caused by disaster, repair or modifications done by anyone other than the company. the company does not warrant that the product will not infringe the intellectual property rights of a third party either alone or in combination with other products or in the operation of any process. The company's total liability for breach of these warranties shall be limited to the gross purchase price for the product as stated in the relevant invoice or order under which the activity giving rise to liability arises. Under no circumstances shall the company be liable for indirect, incidental, consequential, special or punitive damages (including loss of business or profits) under any theory of liability hereunder.

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

LIMITATION OF LIABILITY

Notwithstanding anything to the contrary contained herein, the liability of the company under these terms and conditions (whether by reason of breach of contract, tort, indemnification, or otherwise, but excluding liability of the company for breach of warranty (the sole remedy for which shall be as provided under this section) shall not exceed an amount equal to the lesser of (a) the total purchase price theretofore paid by buyer to the company with respect to the product(s) giving rise to such liability or (b) one hundred thousand dollars (\$100,000), whichever is less.

With respect to any software products incorporated in or forming a part of the Product hereunder Buyer intends and agrees that such software products are being licensed and not sold, and that the words "purchase", "sell" or similar or derivative words are understood and agreed to mean "license", and that the word "Buyer" or similar or derivative words are understood and agreed to mean "licensee". Notwithstanding anything to the contrary contained herein, The COMPANY or its licensor, as the case may be, retains all rights and interest in software products provided hereunder. See the Terms & Conditions of purchase for details regarding such license.

II. Safety



! DANGER (may cause serious damage to property and or casualties)

1. Please carefully read this User Manual prior to operating the instrument and observe the instructions on safe operation.
2. The power source must be grounded reliably and away from any sources of electromagnetic interference.
3. Confirm that the voltage and frequency of the power supply matches the specified voltage and frequency specified prior to use.
4. The instrument should be operated on an independent power outlet.
5. Do not allow the instrument to be plugged in or unplugged whilst the power switch is in the on position.
6. Do not modify the power cord.
7. The working environment must be free of any flammable, explosive, volatile, or corrosive substances.
8. The instrument is de-energized only when the power cable is disconnected from the power source.
9. This instrument is for indoor use only.
10. The instrument should be operated in low humidity, well ventilated, dust and water-free environment, without direct sunshine or an intense magnetic, light, or heat source.
11. Ensure that the instrument and its accessories are free of any potential defects.
12. Ensure only standard accessories supplied by the manufacturers are used with this instrument.

! WARNING (may cause property damage or personal injury)

1. Prevent water from splashing on the electrical elements of this instrument.
2. Before assembly, disassembly, cleaning, or maintenance, the power source must be disconnected.
3. Ensure the vacuum pump's power is not more than the maximum power output of the instrument.

! ATTENTION (may affect operational performance or service life)


1. When pulling the power plug out from a power source, do not pull out using the power cord directly.
2. Only accredited and qualified professional repair technicians can open the instrument or conduct required repairs. Persons performing repairs on the instrument other than those selected or approved by the Company shall operate to void any warranty contained hereinabove for the product

III. Introduction

The instrument is intended to control the vacuum level of a vacuum system with gas as the medium in non-residential environments such as schools, laboratories and or factories.

Environment requirements are as follows:

- Voltage fluctuation does not exceed nominal value by $\pm 10\%$.
- Place the instrument 30cm away from the wall and at least 10 cm from any other instruments.



3.1 After Sales Support

If problems are encountered or technical support is required when installing or using the instrument, please contact ServiceEurope@hollandgreenscience.com.

The Company may provide technical assistance and information regarding the instrument or equipment or service without charge at its sole discretion. Buyer assumes sole responsibility for any reliance on or use of such assistance and information, and the Company makes no warranty thereon.

Upon contact the following information is required:

- Product serial number (located on the instrument nameplate)
- Description of issue or problem
- Method and or operating steps you have undertaken towards resolution.
- Your contact details inclusive of telephone number and email address.

3.2 Proper Use

The instrument is designed for non-residential use and to be used in conjunction only with accessories recommended within this manual and by the manufacturer.

IV. Technical Specifications

Model	20209001
Product Name	Aeolus
Control Range	1-1000mbar
Measurement Range	1-1000mbar
Control Operation Modes	Single-point control operation mode, program operation mode
Control Program	Up to 5 Programs. Each program contains five-step control (vacuum degree and time can be set).
Interface	Touch Screen setting
Display Screen	5-inch TFT
Sensor pressure limit	1500mbar
Connecting Tube Interface Diameter	6.2mm
Materials in Contact with Process	PTFE, PP, Silica Gel and Ceramics
Process Temperature	5°C to 40°C
External Dimension W x D x H	189mm×207mm×193mm
Weight	4.0kg
Input Voltage	AC110—240V, 50Hz/60Hz
Ambient Temperature	10°C-40°C, ≤80%RH
Input Power	600W
Maximum Output (Control) Power	500W(220V), 300W(110V)

V. Inspection

Packing List

Unpack the equipment carefully and check for any damage which may have arisen during transport. In the event of identified damage, please contact ServiceEurope@hollandgreenscience.com.

The package includes the following items:

Item Description.	Quantity.
Main Unit	1
Power Cable	1
Quick Setup Card	1




Note:

If you find any visible damage on the product, please do not connect the product to the power supply.

VI. Operating Instructions

6.1 Overview





Name	Description
Touch Screen	Setting parameters and displaying settings and operation information
Three-in-One Switch	Unit power input, with a switch controlling the power supply to the entire unit.
Pump Connector	Connects the vacuum pump
Air Vent	Connects with air or inert gas if necessary.
Application Connector	Connects with the application.
Power Output Socket	Provides power (control) to vacuum pump.

6.2 Instrument Preparation

- Place the instrument on a stable and well-ventilated surface, clear of flammables or explosives.
- Ensure the instrument voltage matches the power supply voltage before connection.
- Ensure that the power socket is grounded properly.
- Connect vacuum hose of pump and application.
- Connect the power cable of the vacuum pump to the instrument and turn 'ON' the vacuum pump switch. Do not exceed the rated output power of the controller.
- Connect the power cable of the instrument and turn 'ON' the power switch.

6.3 Start-up Interface Description



Button Name	Function
VENT	Press VENT to vent and release to stop venting.
PROGRAM	Choose the control mode and set the control pressure value and running time.
SETTINGS	Set the delay time to 'open and close' the pump and valve. Set the local atmospheric pressure and then adjust the Pressure sensor calibration parameters.
START	Press START to commence the control program, the icon will change to  , press to 'STOP' the control program.
	'PROGRAM' and 'SETTINGS' buttons cannot be operated if the control program starts.


6.5 Operating Mode - Selection and Configuration

Press the 'PROGRAM' button to display the 'Program Selection Interface'




6.6 Single-Point Control Operation Mode


Press 'VACUUM CONTROL' refresh the current page



Name	Function
Set Vacuum (mbar)	Set the vacuum value within a range of 1 - local atmospheric pressure (P0)
Hysteresis (mbar)	Hysteresis value affects the pressure fluctuation range and the 'ON/OFF' frequency of the vacuum pump. Setting range 1-300mbar
Time (min)	Time can be set within a range of 0-5999 min

After setting the parameters, press the 'OK' key to return to the 'Start-Up Interface'.

	<p>Note:</p> <p>The larger the hysteresis value set, the larger the pressure fluctuation, the smaller the hysteresis value set, the smaller the pressure fluctuation and the more frequently the vacuum pump starts up and shuts down. The recommended hysteresis value is set at 10mbar.</p>
---	--

	<p>Note:</p> <p>The sum of the set vacuum value and hysteresis value cannot exceed (P0-10) mbar and the set vacuum value is larger than the hysteresis value.</p>
---	--

6.7 Program Operation Mode

Press 'VACUUM CONTROL' button to display the 'Program Selection Interface'

Function	Program	
Vacuum control	Open :	<input type="text" value="2"/>
Program	Hysteresis(mbar) :	<input type="text" value="10"/>
	<input type="button" value="Edit"/>	<input type="button" value="OK"/>


Name	Function
Open	Select the program number
Hysteresis (mbar)	Hysteresis value affects the pressure fluctuation range and the 'ON' frequency of the vacuum pump Setting range 1-300mbar

After setting the parameters, press the 'OK' key to enter the program operation

6.8 Program Editing


Press the 'Edit' button in the 'Program Selection Interface' to display the 'Program Editing Interface'

Program Editing Area



Name	Function
NO	Sequence No. of program execution
Time (min)	Timing, set within a range of 0-5999min
Vac (mbar)	Set the vacuum value within a range of 1 - local atmospheric pressure (P0)

After setting the parameters, press the 'Program Number' to save until the upper part of the interface indicates 'SAVE'



Note:

If the set vacuum value or timing is blank, it indicates that the current step will be skipped to the next control program.

6.9 Ventilation

Press the 'VENT' button to reduce vacuum.

VII. Fault Diagnosis.

Fault	Cause Analysis	Solution
The vacuum pump does not work.	<p>Check the power switch on the vacuum pump is 'ON'.</p> <p>Check the power cable connection.</p>	<p>Switch the vacuum pump, power switch 'ON'.</p> <p>Connect the cable.</p>
Error message displays when setting the pressure value	Check whether the SET value exceeds the limit.	Adjust the SET value.
The solenoid valve and pump switches ON / OFF frequently.	<p>Check the pipeline for airtightness.</p> <p>Check whether the hysteresis value is too low.</p>	<p>Connect the pipeline properly.</p> <p>Increase the hysteresis value.</p>

VIII. Maintenance and Cleaning

Keep the instrument dry and clean in routine operations.

Clean the outer surface with a non-abrasive cleanser and only connect the power supply when the entire instrument surface is dry.

If liquid or moist solids enters the instrument, please immediately disconnect the power supply and contact supportusa@hollandgreenscience.com

- Surface stains on the instrument should be cleaned only with a clean, soft rag and detergent.
- The power must be disconnected before any maintenance or cleaning.
- Do not clean the instrument with any corrosive cleaning solutions.
- If the instrument is left unused for a long period, switch 'OFF' the power and store it within a dry, clean, level, and stable surface at normal temperature.



Caution!

Before any maintenance or inspection, the power cable **MUST** be removed from the power socket.

Holland Green Science Europe B.V.

Mauritskazerne D41, Office D21

Nieuwe Kazernelaan 2

6711 JC EDE, The Netherlands

Tel: +31643498156

E-mail: SalesEurope@hollandgreenscience.com

E-mail: ServiceEurope@hollandgreenscience.com

Chamber of Commerce number: 82526494

VAT number: NL862505744B01